



Penumbra
3D



Penumbra 


www.penumbrainc.com

Penumbra 3D

A New Way to Remove Clot

Separator™ 3D offers a new thrombus debulking method to the Penumbra System®. Four intraluminal 3D elements secure clot to maximize revascularization with aspiration. Separator 3D is designed to:

- Secure thrombus in the third, or radial, dimension
- Pose minimal vessel wall contact and be atraumatic
- Address the unique challenges of ischaemic stroke



Separator™ 3D*

Catalog Number	Description	Device Diameter (mm)	Device Length (mm)	Wire Length (cm)	Aspiration Catheters
PSS3D	Separator 3D	4.5	26	200	Reperfusion Catheters 041 and 054

* Mikrodrahtretriever German OPS Code 8-83b.8



www.penumbrainc.com

Penumbra, Inc. USA
1351 Harbor Bay Pkwy
Alameda, CA 94502
USA
1.888.BRAIN.06
1.888.272.4606
T 1.510.748.3200
F 1.510.748.3232
order@penumbrainc.com
info@penumbrainc.com

Penumbra Europe GmbH
Am Borsigturm 44
13507 Berlin
Germany
T +49 30 2005 676-0
F +49 30 2005 676-10
order@penumbrainc.de
info@penumbrainc.de

Penumbra Neuro Australia PTY LTD
Level 10, 1 Margaret Street
Sydney NSW 2000
Australia
T +61-1300 817 025
F +61-1300 817 026
order.anz@penumbrainc.com