



Usable Length (80cm)	Usable Length (125cm)	Stent Length (mm)	Stent Diameter (mm)	Sheath
ZX04020080B	ZX04020125B	20	4	6F
ZX04030080B	ZX04030125B	30	4	6F
ZX04040080B	ZX04040125B	40	4	6F
ZX04060080B	ZX04060125B	60	4	6F
ZX04080080B	ZX04080125B	80	4	6F
ZX04100080B	ZX04100125B	100	4	6F
ZX05020080B	ZX05020125B	20	5	6F
ZX05030080B	ZX05030125B	30	5	6F
ZX05040080B	ZX05040125B	40	5	6F
ZX05060080B	ZX05060125B	60	5	6F
ZX05080080B	ZX05080125B	80	5	6F
ZX05100080B	ZX05100125B	100	5	6F
ZX05120080B	ZX05120125B	120	5	6F
ZX05150080B	ZX05150125B	150	5	6F
ZX06020080B	ZX06020125B	20	6	6F
ZX06030080B	ZX06030125B	30	6	6F
ZX06040080B	ZX06040125B	40	6	6F
ZX06060080B	ZX06060125B	60	6	6F
ZX06080080B	ZX06080125B	80	6	6F
ZX06100080B	ZX06100125B	100	6	6F
ZX06120080B	ZX06120125B	120	6	6F
ZX06150080B	ZX06150125B	150	6	6F
ZX07020080B	ZX07020125B	20	7	6F
ZX07030080B	ZX07030125B	30	7	6F
ZX07040080B	ZX07040125B	40	7	6F
ZX07060080B	ZX07060125B	60	7	6F
ZX07080080B	ZX07080125B	80	7	6F
ZX07100080B	ZX07100125B	100	7	6F
ZX07120080B	ZX07120125B	120	7	6F
ZX07150080B	ZX07150125B	150	7	6F
ZX08020080B	ZX08020125B	20	8	6F
ZX08030080B	ZX08030125B	30	8	6F
ZX08040080B	ZX08040125B	40	8	6F
ZX08060080B	ZX08060125B	60	8	6F
ZX08080080B	ZX08080125B	80	8	6F
ZX08100080B	ZX08100125B	100	8	6F
ZX08120080B	ZX08120125B	120	8	6F
ZX08150080B	ZX08150125B	150	8	6F
ZX09020080B	ZX09020125B	20	9	6F
ZX09030080B	ZX09030125B	30	9	6F
ZX09040080B	ZX09040125B	40	9	6F
ZX09060080B	ZX09060125B	60	9	6F
ZX09080080B	ZX09080125B	80	9	6F
ZX09100080B	ZX09100125B	100	9	6F
ZX09120080B	ZX09120125B	120	9	6F
ZX09150080B	ZX09150125B	150	9	6F



ZENFLEX™
Peripheral Stent System



Peripheral Stent System




Superior Stent Performance

The ZENFLEX™ Peripheral Stent is a laser-cutting Nitinol self-expanding stent. Its superior performance is assured by optimal helical stent design combined with Micromesh design.


Micromesh design

- High radial strength and good fracture resistance
- To maximize vessel patency by decreasing tissue prolapse




Unique Helical Design

- Excellent flexibility
- Superior vascular compliance
- High radial support strength



Gold/Tantalum Marker

- 6 embedded micro markers at each end of the stent
- Accurately mark stent ends out and enhance visibility



Controlled Release System

- One-handed stent deployment
- The initial precise positioning followed by well matching quick release
- Ergonomic handle with anti-skid tuning dial combines comfortable hand feel and placement accuracy
- Guidewire acceptance 0.035"



Optimal Flexibility



ZYLOX

The stent is flexible, and both radial force and rate of fracture (<1%) are good, which ensure blood flowing smoothly.



Competitor A

The radial force of the stent is weak and rate of fracture is in the middle level (5-7%). The stent has large lacuna and strong stimuli to vessel wall.



Competitor B

The stent is not flexible with rate of fracture (11-15%), and blood will be hindered in the bending site.

