

CGuard Carotid Embolic Prevention System

SYSTEM SPECIFICATIONS

Size:	Diameter	6 mm-10 mm	✔️ 6 French delivery system
	Length*	20 mm-60 mm	✔️ Self-expanding
Guiding Catheter Compatibility		8F (ID: >2.20 mm or 0.086")	✔️ Rapid exchange
Vascular Sheath Compatibility		6F (ID: >2.20 mm or 0.086")	✔️ Minimal foreshortening
Rapid exchange (RX) Delivery System		6F (OD: 2.03 mm)	✔️ Highly visible under all modalities
Usable Catheter Length		135 cm	✔️ Allows perfusion to ECA
Guidewire Compatibility		0.014"	✔️ Precise placement accuracy
MicroNet® Material		PET	✔️ Optimal endothelialization
Fiber Size		20 µm	✔️ Excellent conformability
Aperture Size**		150 µm - 180 µm	
Stent Material		Nitinol	
Strut Thickness		240 µm ± 12 µm	

TABLE OF SIZES

		Diameter (mm)				
		6	7	8	9	10
Length (mm)	20	CRX0620	CRX0720	CRX0820	CRX0920	CRX1020
	30	CRX0630	CRX0730	CRX0830	CRX0930	CRX1030
	40	CRX0640	CRX0740	CRX0840	CRX0940	CRX1040
	60*	CRX0660		CRX0860		CRX1060

* 60mm sizes are not available in Australia

** Average in Vessel

InspireMD
www.inspiremd.com

Contact Information:

InspireMD Inc.
Information@inspiremd.com
clinical@inspiremd.com
www.inspiremd.com



Not available for sale in the USA

CGUARD™
Carotid Embolic Prevention System




ENGINEERED TO PREVENT STROKES

CONTINUOUS EMBOLIC PREVENTION
WITH MICRONET® SLEEVE

InspireMD


The MicroNet[®] Advantage



Designed to Prevent Distal Embolization and Stroke During and After CAS

CGuard is covered with a uniform PET MicroNet[®] Sleeve for Continuous Embolic Prevention.

MicroNet[®] is a bio-stable mesh woven from a single strand of 20 μm Polyethylene Terephthalate (PET). The MicroNet[®] Sleeve is designed to trap and seal thrombus and plaque against the vessel wall, preventing embolization.




Win-Win Design


CGuard[™] EPS combines the benefits of open cell flexibility with the innovative MicroNet[®] Sleeve—providing maximum surface protection against potential emboli.

MicroNet[®] Technology


- ✓ Traps and seals emboli
- ✓ Lends greater confidence in post-dilation
- ✓ Prevents plaque prolapse and late embolization
- ✓ No restenosis concerns at 12 months*




CGuard OCT



CGuard 3D OCT



CGuard Radiology



CGuard IVUS




Open Cell Nitinol Structure


SmartFit[™] Technology

Eliminates the need for tapered design.

SmartFit Technology maintains optimal placement length regardless of vessel diameter.



CGuard MicroNet[®] Sleeve has the smallest pore size of any carotid stent for continuous embolic prevention.



*Data on File